



8 Day Kitchen Clock Movement



Movement size 5-1/2" X 6" X 2-5/8" deep
Over all with 1" handshaft. Pendulum drop
Is 9-1/4" from center of handshaft to bottom
of pendulum tip. This movement mounts
to the back of the case. Accessories not
included.

24.810 \$ 52.75 each.


8 Day Mantle Clock Movement



Movement size 5-1/2" X 6" X 3-3/4" deep
Over all with 1-1/2" handshaft. Pendulum
drop Is approximately 4-1/2" from center of
handshaft to bottom of pendulum tip.
This movement is front mount with the
Pendulum hanging in the rear.
Accessories not included.

24.811 \$ 52.75 each.

Gathering Pallets



79.342 Left



79.343 Right

Used on Tall Case clocks.


79.342 \$ 18.50 ea.

79.343 \$ 18.50 ea.

Mix or match 2 or more ... \$ 18.00 ea.

Prices Subject To Change Without Notice.

8 Day School Movement



This movement can be used in School
Clocks, long drop store regulators or
regulators clocks.

Movement size 5-1/2" X 6" X 3-3/4" deep
Over all with 2" hand shaft length.

Pendulum drops are measured from center
of hand shaft to bottom of pendulum tip.
Accessories not included.

24.812 ... 13-1/4" drop. ... \$ 52.75 each

24.813 ... 18" drop \$ 52.75 each.

Movement Test Stand



Hold clock movements when regulating
and making adjustments. Grooved brackets
are adjustable up to 11". 3 leveling screws
are in the base. Post diameter 1/2".

24.800 \$ 39.50 each

Movement Test Stand

Aluminum test stand
Completely adjustable.

It has three thumb turns
For leveling in all direction

The horizontal arms are


grooved to hold
movements. The

round collors are
designed to hold
round movements.


The rugged

aluminum
stand. The stand is 16-1/2" tall with a 1/2"
post diameter.

24.801 \$ 43.50 ea.



Alarm Mechanism



Alarm mechanism with center alarm ring
And actuating arm. For Kitchen Clock and
Etc. Comes with mechanism and hammer.

24.809 \$ 12.50 ea.



Center Alarm ring with lever.

24.808 \$ 4.00 ea.


MOVEMENT HOLDER



Movement holder can be used to assemble
movements and then hang on a board or the
wall for testing. Included in the set is 4 post
1 Strap Bracket And a wall stud for the
holder to hang on.

28.802 \$ 15.50

Long Movement Assembly Post



Assembly Post 3 3/4" long. Set of 3

28.805 \$ 16.25 set

28.806 ... set of 4 pc. ... \$ 21.75 set